# COVID-19 VACCINES: TEMPERATURE MONITORING PLUG&TRA



#### 2 SOLUTIONS FOR VACCINE STORAGE TEMPERATURE CONTROL

## DESCRIPTION

You are receiving Covid-19 vaccines that you need to store in your refrigerators. Once the Pfizer/BioNTech vaccine is delivered to the vaccination point the storage time is 5 days at a temperature between +2 and +8°C. The Moderna vaccine should be stored at -20°C, the others between +2 and +8°C.

In order to comply with the HAS guidelines and to ensure the temperature control of the vaccines, we propose two solutions:

- A traceability system with a temperature data logger to be read periodically.
- An automated alarm and traceability system in real time.



### **SOLUTION 1: TRACEABILITY SYSTEM**

- 1 Temperature data logger Thermo Button: -40/+85°C ±0.5°C
- 1 Thermotrack PC Solo software and reader
- Recommended option : calibration certificate

Requires periodic reading on a PC

125€

200€






Software and Thermo Button reader

## **SOLUTION 2: AUTOMATED ALARM** AND TRACEABILITY SYSTEM

- 1 Digital temperature sensor : -40 à +85°C ±0.5°C
- 1 Sensor Net Connect Ethernet smart hub
- 1 Thermotrack Webserve Annual Alarm and **Traceability System**

**Automatic real-time readings** 



BOTH SOLUTIONS ARE COMPATIBLE WITH ALL VACCINES: PFIZER, MODERNA AND OTHERS

#### **EQUIPMENT REQUIRED ACCORDING TO THE CDC:**

Storage unit temperatures must be monitored regularly and checked and recorded at the beginning of each workday to determine if any excursions have occurred since the last temperature check. For accurate temperature monitoring, use a digital data logger (DDL)"

Plug and Track A Proges Plus activity Lille - Lyon - Nantes 2, rue de la République 59780 - WILLEMS - France Phone: + 33 3 20 64 63 63 info@plugandtrack.com www.plugandtrack.com

Pfizer-BioNTech COVID-19 Vaccine Storage and Handling Summary Centers for Disease Control and Prevention, USA, december 2020



#### I CHOOSE MY SOLUTION:

/	REFERENCE	DESCRIPTION	QUANTITY	UNIT PRICE €	TOTAL	
	TRACEABILITY SYSTEM					
	PAC0044	Pack Solo 22L with Thermo Button 22L, Thermotrack PC Solo software and reader : -40/+85°C ±0.5°C		182.00€		
	ZA001	Option : calibration certificate		50.00€		
		AUTOMATED ALARM AND TRACEAR	BILITY SYSTE	M		
	PAC0089	Pack Pharma Fridge with Sensor Net Connect Ethernet temperature sensor and annual subscription to Thermotrack Webserve		278.00€		
	ZA0001	Option : calibration certificate		50.00€		
					,	
	PR20017	Shipping costs : France 20€/Europe 40€				
	For the European	For the European Union, please indicate your VAT number				
	YOUR COMPANY					
			Name :			
	Company stamp	:	Date:			
			Signature :			

**Delivery time:** 24 hours without calibration. 2 to 3 weeks with calibration.

TERMS AND CONDITIONS OF PAYMENT: to the order by check or bank transfer. IBAN: FR76 3002 7175 2800 0241 8060 129

**DELIVERY NOTE**: to be sent by email to thermotrack@proges.com